



# YOUR INDIVIDUAL PRODUCTS CATALOG

Creation date 21.09.2024





Creation date 21.09.2024, Subject to change without notice


[Click here](#)

[To the shop](#)



## POSITIONING PLATE SET

Item No: H06-561.002



### TECHNICAL DATA



### PRICE

#### On request

All prices net, plus VAT and delivery  
Valid on 21.09.2024, Changes and errors reserved

### SPECIFICATIONS

- **upholstery**  
black, antistatic
- **cushion cover**  
hard-wearing, protective coated artificial leather
- **cushion cover**  
hard-wearing, protective coated artificial leather

### PRODUCT DESCRIPTION

The support plates suitable for the MultiPro® 2 / 5 treatment table can be aligned on both sides. The set consists of a pair of reinforced equipment rails, a back section pad, a leg section pad and an extendable leg section pad. The assembly requires a pad width of 65 cm and a pad thickness of 5 cm. With this positioning aid, the patient can be positioned easily and safely in order to perform a rectal examination.

