



YOUR INDIVIDUAL PRODUCTS CATALOG

Creation date 21.09.2024






Creation date 21.09.2024, Subject to change without notice

[Click here](#)

[To the shop](#)



Item No: H20-189.01



TECHNICAL DATA

illuminance
50'000 Lux

continuous dimming
20 - 100 %

light field Ø
26 cm

lamp body Ø
19.5 cm

colour rendering index
95 RA

colour temperature
4000, 4500 Kelvin

illuminance
50'000 Lux

continuous dimming
20 - 100 %

light field Ø
26 cm

lamp body Ø
19.5 cm

colour rendering index
95 RA

colour temperature
4000, 4500 Kelvin

PRICE

On request

All prices net, plus VAT and delivery
Valid on 21.09.2024, Changes and errors reserved

PRODUCT DESCRIPTION

The HESS MultiLED 9 led operating lamp, which generates no shadows, is suitable for small surface surgeries, gynaecology and first aid. The close arrangement of the optics (the reflector size is 195mm) does not require focusing. The lamp is easy to handle thanks to the aluminum supporting structure, which gives the product lightness in movements. The reduction of the LED's IR rays, together with the low power used, allow to obtain a low heat emission, for the benefit of the patient but

SPECIFICATIONS

- **execution**
ceilings mounted
- **intensity adjustment**
by means of a control panel on the front of the lighting body
- **positioning**
by means of a central, sterilizable handle
- **Configuration**
ceiling mounting
- **frame colour**
signal white RAL 9003
- **surface**
easy to clean & disinfect



also of the doctor.





Creation date 21.09.2024, Subject to change without notice

[Click here](#)

[To the shop](#)



Item No: H20-189.01

VARIANTS

You will find further variants in our Online Shop



**Untersuchungsleuchte HESS
MultiLED 12 Decke**

H20-188.01

