



YOUR INDIVIDUAL PRODUCTS CATALOG

Creation date 21.09.2024





Creation date 21.09.2024, Subject to change without notice

[Click here](#)

[To the shop](#)



THERAPY CHAIR WITH SADDLE SEAT

Comfort for your everyday life

Item No: H01-153.009



TECHNICAL DATA

cushion thickness
6 cm

height adjustment range
48 - 66 cm

double swivel castors
Ø 50 mm

maximum load capacity
150 kg

amount of spokes
5 Stk.

foot cross
Ø 57 cm

amount of spokes
5 Stk.

foot cross
Ø 57 cm

PRICE

On request

All prices net, plus VAT and delivery
Valid on 21.09.2024, Changes and errors reserved

COLOURS



PRODUCT DESCRIPTION

The high-quality therapy chair not only offers you comfort. The seat, covered with high-quality artificial leather, is extremely durable and resistant to disinfectants. The medical upholstery not only offers a pleasant sitting sensation, it is also seamlessly welded air- and watertight. No compromises when it comes to hygiene in your medical practice!

SPECIFICATIONS

- **upholstery**
HESS «MediPad» cushion: protective coated and fully sealed upholstery, therefore absolutely leak-proof against all liquids
- **castor type**
plastic, stainless with thread guard
- **Seat area**
saddle seat
- **height adjustment**
by means of hand release lever
- **foot construction**
pyramid-shaped 5 star foot cross, stainless steel ground and electropolished
- **SWISS MADE**
a Swiss product manufactured by HESS
- **Configuration**
48 - 66 cm, double swivel castors





Creation date 21.09.2024, Subject to
change without notice

[Click here](#)

[To the shop](#)



THERAPY CHAIR WITH SADDLE SEAT

Comfort for your everyday life

Item No: H01-153.009

ACCESSORIES

You can find further accessories in our Online Shop



HESS MultiClean®
detergent 5 liter

H06-991.01



HESS MultiClean®
detergent 300 ml

H06-991

