



YOUR INDIVIDUAL PRODUCTS CATALOG

Creation date 21.09.2024





Creation date 21.09.2024, Subject to change without notice

[Click here](#)


[To the shop](#)



HESS SOFT STRETCH MATTRESS

Suitable for the HESS C-arm and X-ray table SwissTroll® X

Item No: H06-558.014



TECHNICAL DATA

Outer dimensions L x W x D
ca. 210 x 77 x 6 mm

patient weight
60-100 kg

PRICE

On request

All prices net, plus VAT and delivery
Valid on 21.09.2024, Changes and errors reserved

PRODUCT DESCRIPTION

The HESS mattress is perfectly matched to the HESS C-arm and X-ray table SwissTroll® X. Thanks to the innovative single-core mattress construction, the mattress offers excellent lying comfort. The breathable mattress cover can be changed easily and hygienically and not only guarantees freshness, but is also anti-allergic and absolutely impermeable to body fluids. The mattress is also washable, ensuring optimum hygiene. The mattress is ideal for use in emergency departments and offers the perfect combination of comfort, ease of care and reliability.

SPECIFICATIONS

- **SWISS MADE**
a Swiss product manufactured by HESS
- **Mattress core**
comfortable double core
- **Mattress cover**
imitation leather electrically conductive, bacteria and virus-proof, breathable, impermeable to liquids
- **Cover underside**
underside made of anti-slip material, coated PES fabric with sewn-on fleece tape as optimum anti-slip protection for a secure hold
- **Replacement cover**
the mattress cover can be easily removed by means of a concealed zip fastener removed from the mattress core, or replaced
- **Color of cover**
night blue
- **care**
cover washable up to 95° C, drying at low temperature in tumbler possible
- **Disinfection**
wipe disinfectable with disinfectant according to recommendation list

