



# YOUR INDIVIDUAL PRODUCTS CATALOG

Creation date 21.09.2024





Creation date 21.09.2024, Subject to change without notice

[Click here](#)

[To the shop](#)



# THERAPY CHAIR BASELINE WITH HAND RELEASE

Comfort for your everyday life

Item No: H01-320.300



## TECHNICAL DATA

upholstery dimensions  
Ø 36 cm

cushion thickness  
6 cm

height adjustment range  
51 - 69 cm

double swivel castors  
Ø 50 mm

maximum load capacity  
150 kg

amount of spokes  
5 Stk.

foot cross  
Ø 50 cm

amount of spokes  
5 Stk.

foot cross  
Ø 50 cm

## PRICE

CHF 389.00 / Stk

All prices net, plus VAT and delivery  
Valid on 21.09.2024, Changes and errors reserved

## COLOURS



## PRODUCT DESCRIPTION

The high-quality therapy chair not only offers you comfort. The seat, covered with high-quality artificial leather, is extremely durable and resistant to disinfectants. The medical upholstery not only offers a pleasant sitting sensation, it is also seamlessly welded air- and watertight. No compromises when it comes to hygiene in your medical practice!

## SPECIFICATIONS

- **upholstery**  
HESS «MediPad» cushion: protective coated and fully sealed upholstery, therefore absolutely leak-proof against all liquids
- **castor type**  
plastic, stainless with thread guards safety castors, locked without load
- **Seat area**  
round seat
- **height adjustment**  
via hand release ring
- **foot construction**  
polished aluminium base frame
- **Delivery**  
pre-assembled
- **Configuration**  
51 - 69 cm





Creation date 21.09.2024, Subject to change without notice

[Click here](#)

[To the shop](#)



## THERAPY CHAIR BASELINE WITH HAND RELEASE

Comfort for your everyday life

Item No: H01-320.300

### ACCESSORIES

You can find further accessories in our Online Shop



HESS MultiClean®  
detergent 5 liter

H06-991.01



HESS MultiClean®  
detergent 300 ml

H06-991

