



# YOUR INDIVIDUAL PRODUCTS CATALOG

Creation date 21.09.2024





Creation date 21.09.2024, Subject to change without notice

[Click here](#)

[To the shop](#)



# PAIR OF SIDE RAILS FOR PATIENT BEDS WITH EQUIPMENT RAIL

Safety for your patients!

Item No: H06-544



## TECHNICAL DATA

length  
110 cm

height  
32 cm

equipment rail  
25 x 10 x 1000 mm

equipment rail  
25 x 10 x 1000 mm

## PRICE

On request

All prices net, plus VAT and delivery

Valid on 21.09.2024, Changes and errors reserved

## PRODUCT DESCRIPTION

The side rails for our patient tables create additional safety for your patients. Made of stainless steel with a robust and high-quality release mechanism – easy and safe to operate without any clamping points. The side rails can be lowered to below the level of the upholstery for optimal access to the table. The integrated equipment rail allows you to attach additional accessories to the table.

## SPECIFICATIONS

- **operation**  
by gently pulling the safety lock, automatically locking upwards
- **execution**  
1 pair, left and right version
- **material**  
made entirely of stainless steel, ground and electropolished
- **mounting type**  
firmly screwed
- **compatibility**  
for upholstery width 65 - 80 cm
- **SWISS MADE**  
a Swiss product manufactured by HESS
- **Configuration**  
with equipment rail

