



YOUR INDIVIDUAL PRODUCTS CATALOG

Creation date 21.09.2024





Creation date 21.09.2024, Subject to change without notice

[Click here](#)

[To the shop](#)



HEIGHT-ADJUSTABLE OXYGEN CYLINDER HOLDER

For 2-, 3- or 5-litre bottles

Item No: PR-Z-DOCK80



TECHNICAL DATA

maximum load capacity
12 kg

capacity
2 - 5 l

PRICE

On request

All prices net, plus VAT and delivery
Valid on 21.09.2024, Changes and errors reserved

PRODUCT DESCRIPTION

Our oxygen supply solution offers decisive advantages for the care and transfer of patients. The continuous oxygen supply from the oxygen cylinder is not interrupted in an emergency or during transfers, ensuring a constant and reliable supply. The technology is characterised by a high-quality coating that ensures a long service life. There are no sharp edges, which minimises the risk of injury and increases safety for patients and nursing staff. The application is particularly user-friendly, as the solution can be retrofitted and is easy to install. This means it can be quickly and easily integrated into existing docking trolleys. Experience a safe, efficient and user-friendly oxygen supply solution.

SPECIFICATIONS

- **Mounting**
Retaining claws with quick-release lever
- **safety**
DIN EN 60601
- **mounting type**
Can be fixed on equipment rails 25 x 10 mm
- **care**
easy to clean by wipe disinfection

